



**WARDRAY
PREMISE**

Economy Fixed Height Trolley

MR1510

Tested for use in MR environment up to 7 Tesla

Table top dimensions 1930mm x 550mm

Overall dimensions including cot sides 2100mm x 625mm x 800mm

Backrest adjustable to 85°

Cot sides can be lowered

Fitted with 150mm braked swivel castors

Powder coated in standard MRI blue

50mm comfort mattress in grey with stitched seam

Weight capacity up to 220Kg (490lbs)



Optional extras:

IV Pole

Patient restraining steps

Oxygen Cylinder Holder

Cot bumpers

GEN.S.16.295