



FOLDING BARIATRIC PORTERING CHAIR: MR4588

MR



Dimensions:

Overall:

Height: 925mm
Width: 742mm
Depth: 1055mm
Width (folded): 475mm

Seat:

Height: 550mm
Width: 630mm
Depth: 550mm

This new chair has been designed specifically to transport bariatric patients to and from the MRI suite and is MR safe.

With the increase of obesity becoming more of an issue in recent years, this chair allows patients to be transported safely to an open scanner.

The MR4588 is approximately 30% wider than our standard portering chair and has been weight tested to 300kg (approx 47 stone) static weight.

With no magnetic parts included in the design, the chair can transport the patient in to the scanner room.

Armrests raise and lower to below seat level to ease patient transfer.

**THE MR4588 IS FULLY
TESTED MR SAFE**

Wardray Premise Ltd

Hampton Court Estate
Summer Road
Thames Ditton
Surrey
Tel: +44 (0) 20 8398 9911
Fax: +44 (0) 20 8398 8032
E-mail: sales@wardray-premise.com

Full details of all products and services can be viewed on our website:

www.wardray-premise.com

QUALITY WITHOUT
COMPROMISE SINCE 1909